

The Biocon exosome separation platform is a method and instrument to separate exosomes based on size difference without any kind of force which can damage the exosome. Undamaged exosomes are important to study intercellular cell signaling process which are based on functionality of exosomes. We suggest two kinds of separation model, one is based on Brownian motion and the other is using electrophoretic mobility difference of particle. The Biocon exosome separation platform can play a key role in investigate biochemical cell signaling process.

(Copyright © 2012, American Society for Microbiology. All Rights Reserved.)

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
(Copyright © 2012, American Society for Microbiology. All Rights Reserved.)

### Exosome Separation by Brownian Motion

- Chip**
- Chip is composed of PDMS and PC (polycarbonate) membrane filter
  - This chip can separate exosomes by size
  - This technique is based on Brownian motion, it does not give damage to exosomes after separation
- Separation**
- Between the two PDMS chip, nanoporous membrane filter is placed which has specific pore size (50nm, 100nm, 200nm, 400nm)
  - Exosomes which are smaller than membrane pore size can pass through the filter

### Exosome Separation by Lateral Electrophoresis


- Chip**
- Material : B270 glass
  - AZ 1500 PR for etching mask
  - Specific pattern can be made by film mask and own lithography device
  - Wet etching using 10:1 BOE
  - Chip size is 1X1 inch and height is 1mm
  - Channel width and height are 60, 12 um

- Separation**
- Electrophoretic mobility difference based separation which is determined by particle size and surface charge
  - Electric field is applied to sample flow(main channel) by side channel which has a role to supply ions to main channel through ion selective polymer
  - Both channels are separated by ion selective polymer which exchange only ions

### ARS enzyme kinetics chip

- Chip**
- PDMS chip
  - It is designed to measure sample quantity, make concentration gradient and mix 36 kinds of parallel reaction can be performed simultaneously
  - Reagents induced to channel can be quantized by the distance of horizontal valve
  - Each processors are designed to mix different concentration ratio
  - Fast fluid mixing using peristaltic mixing valves

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